



# PRODUCT SPECIFICATION



**Product:** CEDERROTH 4-IN-1 MINI

**Reference code:** 1911



## 1. General description

Cederroth 4-in-1 Mini is a sterile universal dressing, which consists of one elastic bandage and one compress with a non-stick wound contact layer which is stitched to the bandage. The dressing is wrapped in a pouch with printed pictorial user instructions.

## 2. Intended use

To be used in first aid situations. Cederroth 4-in-1 Mini is especially suitable for fingers and toes has four functions:

- Pressure pad for heavy bleeding
- Protective dressing for superficial wounds
- Protective dressing for superficial wounds
- Burn dressing with special wound protecting surface

Because of the many functions, Cederroth 4-in-1 replaces a variety of traditional first aid dressings and compresses.

## 3. Technical description

### *Product dimensions :*

Consumer package: 6,5 x 4 x 2,5 cm  
Bandage: 300 x 6 cm  
Compress: 12 x 8 cm

### *Material*

Bandages Rayon  
Compress Absorbent inner layer: Rayon staple fibers  
Non adherent outer cover: Porous net made of high density polyethylene  
Packing material (sterile package) Transparent film made of polyethylene and polyamide. Paper.

## 4. Sterilization

Sterilization method: Irradiation

## 5. Composition

### *Consumer package*

### *Transport package*

<u>Ref#</u>	<u>Product/Material</u>	<u>Pcs</u>	Number of consumer packages	
51616560	Small bloodstopper	1	in transport package:	30
51616575	Plastic pouch	1	Dimensions (l x w x h) cm:	21 x 16 x 11
			Weight (kg):	0,5



## PRODUCT SPECIFICATION



### 6. Labeling

*On product/ consumer package*

Cederroth address, telephone- and faxnumber, web-address, mail-address

Product name, Reference code

Lot-number, Expiry date, EAN-code,

Single use, Sterile R, User

instructions, Product size, CE+NB#

*On transport package*

Cederroth address, telephone- and faxnumber, web-address

Product name, Reference code

Lot-number, Expiry date, EAN-code, Quantity, CE+NB#

### 7. Storage conditions

The product should be stored in the original package for maximum protection until time of use. Do not use if the sterile package is broken.

Shelf life: 5 years from production date

### Authorization: