



# PRODUCT SPECIFICATION



**Product:** SALVEQUICK WOUND CLEANSER  
REFILL  
**Reference code:** 323700



## 1. General description

Sterile sachets made of non-woven tissues, which are saturated in a wound cleansing liquid. The tissues are packed in a well-sealed foil, which prevents leakage and guarantees sterility. The sachets are packed in a refill box.

## 2. Intended use

Salvequick Wound Cleanser is intended for mechanical cleansing of minor wounds when clean tap water is not available.

## 3. Technical description

### *Product dimensions*

Sachet	(l x w ) cm:	6 x 6
Tissue	(l x w ) cm:	12 * 11
Box	(l x w x h) cm:	6 x 6 x 8

### *Material*

Tissue:	White non-woven tissue
Foil:	Laminate of Paper/Aluminum/Surlyn® with print
Liquid:	Physiological saline solution

## 4. Sterilization

Sterilization method: Irradiation

## 5. Composition

### *Consumer package*

<u>Ref#</u>	<u>Product/Material</u>	<u>Pcs</u>
3239	Salvequick Wound Cleanser sachets	20

### *Transport package*

Number of consumer packages in transport package:	20
Dimensions (l x w x h) cm:	33 x 14 x 16
Weight (kg):	1,9



# PRODUCT SPECIFICATION



## 6. Labeling

### *On product (sachet)*

Cederroth address, telephone- and faxnumber, web-address  
Product name, Reference code  
Lot-number=expiry date, Expiry date  
Symbols: CE +NB#, STERILE R, DO NOT REUSE, LOT, USE BY

### *On consumer package*

Cederroth address, telephone- and faxnumber, web-address  
Product name, Reference code  
Lot-number, Expiry date, EAN-code  
Symbols: CE +NB#, STERILE, LOT, USE BY

### *On transport package*

Cederroth address, telephone- and faxnumber, web-address  
Product name, Reference code  
Lot-number, Expiry date, EAN-code  
Quantity  
Symbols: LOT, USE BY

## 7. Storage conditions

The product should be stored in the original package for maximum protection until time of use.

Shelf life: 5 years from production date

## 8. Accessories

No accessories are needed for the use of this product.

## Authorization: